

重组蛋白

Catalog Number:ORSEP5404AMu

Recombinant Escherichia coli Iron-sulfur cluster repair
protein ytfE (ytfE) (Active)

Organism Species: Escherichia coli (strain K12)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: YTFE

Uniprot:P69506

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: Not test.

Source: E.coli

Protein Length: Full Length

Expression Region: 1-220aa

Sequence:

MAYRDQPLGELALSIPRASALFRKYDMDYCCGGKQTLARAAARKELDVEVIEAELAKLAE
QPIEKDWRSAPLAEIIDHIIIVRYHDRHREQLPELILQATKVERVHADKPSVPKGLTKYLT
MLHEELSSHMMKEEQILFPMIKQGMGSQAMGPISVMESEHDEAGELLEVIKHTTNNVTPP
PEACTTWKAMYNGINELIDDLMDHISLENNVLFPRALAGE

Tag Info: N-terminal GST-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.

If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:51.9 kD

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



lyophilized form is 12 months at -20°C / -80°C . Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>

手机号: 13061929181(微信)
13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

