

# 重组蛋白

Catalog Number:ORSEP5407AMu12

Recombinant Human Novel Coronavirus

Nucleoprotein(N) (R203K, G204R) (Active)

Organism Species: Human Novel Coronavirus (SARS-CoV-2/  
2019-nCoV)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

**【Product Information】**

Gene Names: N

Uniprot: P0DTC9

Purity: Greater than 96% as determined by SDS-PAGE.

Endotoxin: Not test.

Source: E. coli

Protein Length: Full Length

Expression Region: 1-419aa (R203K, G204R)

Sequence:

MSDNGPQNQRNAPRITFGGPS DSTGSNQNGERSGARSKQRRPQGLPNNTASWFTALTQHG  
KEDLKFPRGQGVPINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSRWYFYLLGTGPEAG  
LPYGANKDGI IWVATEGALNTPKDHIGTRNPANNAIIVLQLPQGTTLPKGFYAEGSRGGS  
QASSRSSRSRNSSRNSTPGSSKRTSPARMAGNGGDAALALLLDRLNQLESKMSGKGQQ  
QQGQTVTKKSAAEASKKPRQKRTATKAYNVTQAFGRRGPEQTQGNFGDQELIRQGTDYKH  
WPQIAQFAPSASAFFGMSRIGMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLNKHIDAY  
KTFPPTEPKDKKKKKADETQALPQRQKKQQTVTLLPAADLDDFSKQLQQSMSSADSTQA

Tag Info: N-terminal 6xHis-tagged

Form: Lyophilized powder

Buffer formulation: Lyophilized from 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0

Applications: -

Predicted Molecular Mass: -

Accurate Molecular Mass: 49.8 kDa

Predicted isoelectric point: -

**【Storage and Stability】**

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991




企业微信

企业微信

lyophilized form is 12 months at  $-20^{\circ}\text{C}$ / $-80^{\circ}\text{C}$ . Avoid repeated freeze/thaw cycles.

### 【Identification】



### 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com>  
<https://www.aoruicell.com>

手机号: 13061929181(微信)  
13061929770(微信)

邮箱: sale4@aoruisai.com  
sale5@aoruisai.com

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

