

重组蛋白

Catalog Number:ORSMP5411AMu

Recombinant Paguma larvata Angiotensin-converting

enzyme 2 (ACE2), partial (Active)

Organism Species: Paguma larvata (Masked palm civet)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: ACE2

Uniprot:Q56NL1

Purity: Greater than 94.8% as determined by SDS-PAGE.

Endotoxin: Less than 1.0 EU/ug as determined by LAL method.

Source: Mammalian cell

Protein Length: Partial

Expression Region: 18-740aa

Sequence:

QSTTEELAKTFLETFNIEAQELSYQSSVASWNYNTNITDENAKNMNEAGAKWSAYYEEQS
KLAQTYPLAEIQDAKIKRQLQALQQSGSSVLSADKSQRLNTILNAMSTIYSTGKACNPNN
PQECLLLEPGLDNIMENSKDYNERLWAWEGWRAEVGKQLRPLYEEYVALKNEMARANNYE
DYGDYWRGDYEEEWTTGGYNYSRNQLIQDVEDTFEQIKPLYQHLHAYVRACLMDTYPSTRIS
RTGCLPAHLLGDMWGRFWTNLYPLTVPFQKPNIDVTDAMVNNQNDARRIFKEAEKFFVS
VGLPNMTQGFWENSMLTEPGDGRKVVCHPTAWDLGKGDFRIKMCCTKVTMDDFLTAHHEMG
HIQYDMAYAAQPFLLRNGANEGFHEAVGEIMSLAATPNHLKTIGLLSPAFSEDNETEIN
FLLKQALTIVGTLPTFTYMLEKWRWMVFKGAIPKEQWMQKWWEMKRNIVGVVEPVPHDETY
CDPASLFHVANDYSFIRYYTRTIYQFQFQEQALCQIAKHEGPLHKCDISNSTEAGKKLLEM
LSLGRSEPWTALERVVGAKNMNVTPLLNYFEPLFTWLKEQNRNSFVGWDTDWRPYSQS
IKVRISLKSALGEKAYEWNENMYLFRSSIAAYAMREYFSKVKNQTIPFVEDNVWVSDLKP
RISFNFFVTFSNNVSDVIPRSEVEDAIRMSRSRINDAFRLDDNSLEFLGIEPTLSPPYRP
PVT

Tag Info: C-terminal hFc-tagged

Form:Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:112.7

Predicted isoelectric point:-

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991




企业微信

企业微信

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信