

重组蛋白

Catalog Number:ORSMP5413AMu

Recombinant Macaca fascicularis Delta-like protein

3 (DLL3), partial (Active)

Organism Species: Macaca fascicularis (Crab-eating
macaque) (Cynomolgus monkey)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: DLL3

Uniprot:A0A2K5WSR4

Purity: Greater than 82% as determined by SDS-PAGE.

Endotoxin: Less than 1.0 EU/ug as determined by LAL method.

Source: Mammalian cell

Protein Length: Partial

Expression Region: 27-490aa

Sequence:

AGVFELQIHSFGPGPGAPRSPCSARGPCRLFFRVCLKPGLSEEAESPICALGAALSAR
GPVYTEQPEAPAPDLPLPNGLLQVPFRDAWPGTFSLIETWREELGDQIGGPAWSLLARV
TRRRRLAAGGPWARDIQRAGAWELRFSYRARCELPVGTACTRLCRPRSAPSRCGPGLRP
CAPLEDECEAPPVCRAGCSLEHGFCEQPGECCRCLEGWTGPLCMVPVSTSSCLGLRGPSST
TTGCLVPGPGPCDGNPCANGGSCSETPGSFECTCPRGFYGLRCEVSGVTCADGPCFNGGL
CVGGADPDSAYICHCPPGFQGSNCEKRVDRCSLQPCRNGGLCLDLGHALRCRCRAGFAGP
RCEHDLDDCAGRACANGGTCVEGGGAHRCSCALGFGGRNCRERADPCAARPCAHHGRCYA
HFSGLVCACAPGYMGARCEFPVHPDGVSALPAAPPGLRPGDPQR

Tag Info: C-terminal 6xHis-tagged

Form:Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:50.6 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991




企业微信

企业微信

lyophilized form is 12 months at -20°C / -80°C . Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
 sale5@aoruisai.com

电话: 021-57896695
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aorucell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信