

重组蛋白

Catalog Number:ORSYP0323R0

Recombinant *Puccinia triticina* Uncharacterized
protein(PTTG_03405)

Organism Species: *Puccinia triticina* (isolate 1-1 /
race 1 (BBBD)) (Brown leaf rust fungus)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: PTTG_03405

Uniprot:A0A180GV39

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: Yeast

Protein Length: Full Length of Mature Protein

Expression Region: 19-189aa

Sequence:

DAPASSTAATNAKTASAARAAAPGGKQPNLKPMNCTQAFLPFSEADVAAALSSNDTKAAAS
GNYSTIPTEAACKGPAGSIDGLCDIGSCGNHPVCNTCVELIITGPNTTKTGTTTTTPQVTC
TNNYFFGNSTDPIKNVCTDGNDKTYTCTGGCVSYTSCQMCVSVDDPALQGP

Tag Info: C-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:18.7 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552


发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>

手机号: 13061929181(微信)
13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信