

重组蛋白

Catalog Number:ORSYP4837R0

Recombinant *Fusarium solani* subsp. *cucurbitae*

Cutinase (CUTA)

Organism Species: *Fusarium solani* subsp. *cucurbitae*
(*Neocosmosporum cucurbitae*)



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: CUTA

Uniprot:Q99174

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: Yeast

Protein Length: Full Length of Mature Protein

Expression Region: 17-230aa

Sequence:

LPTSNPAQELEARQLGRTTRDDLINGNSASCADVIFIIYARGSTETGNLGTLGPSIASNLE
SAFGTDGVWVIQGVGGAYRATLGDNALPRGTSSAAIREMLGLFQQANTKCPDATLIAGGYS
QGAALAAASIEDLDSAIRDKIAGTVLFGYTKNLQNRGRIPNYPADRTKVFCNVGDLVCTG
SLIVAAPHLAYGPDARGPAPEFLIEKVRVAVRGS

Tag Info: C-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:23.7 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)


邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

