

CD87

产品信息

产品名称	产品货号	产品规格
CD87	ORS7256PAb	50 μ l/100 μ l/1ml

产品简介

CD87, known as urokinase plasminogen activator-receptor (uPAR), is a 36-68 kD single-chain glycoposphotidylinositol (GPI)-anchored glycoprotein. The extracellular region of CD87 consists of three homologous domains, D1, D2, and D3. Membrane-bound and soluble (D2D3) forms of uPAR have been reported. The membrane-bound CD87 is expressed on granulocytes, monocytes/macrophages, dendritic cells, fibroblasts, endothelial cells and keratinocytes. CD87 plays a role in cell adhesion and chemotaxis through binding to β 1, β 2 and β 3 integrins. CD87-bound uPA converts plasminogen to plasmin.

产品性质

Alias: PLAUR, u-PAR, uPAR

Clonality: Recombinant Mouse Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: P35456

Clone: 7A9

Entrez GeneID: 18793

Host\Isotype: Mouse IgG2a

Reactivity: Mouse

Purification: affinity purified by Protein G

Applications: FCM=1/100-1/200

Immunogen: -

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信