

5-HT-6

产品信息

产品名称	产品货号	产品规格
5-HT-6	ORS13887MAb	100 μ l/50 μ l

产品简介

This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of this receptor is mediated by G proteins that stimulate adenylate cyclase. It has a high affinity for tricyclic psychotropic drugs. Kohen R., J. Neurochem. 66:47-56 (1996). Kopatz S.A., Submitted (OCT-2003) to the EMBL/GenBank/DBJ databases. Gregory S.G., Nature 441:315-321 (2006).

产品性质

Alias: 5-HT6; 5HT6R; HTR6;

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: P50406

Clone: -

Entrez GeneID: 3362;

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/500-1/3000

Immunogen: Synthesized peptide derived from internal of human 5-HT-6.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信