

PABPC1

产品信息

产品名称	产品货号	产品规格
PABPC1	ORS8571MAb	100 μ l/50 μ l

产品简介

PABPC1 binds the poly(A) tail of mRNA. It may be involved in cytoplasmic regulatory processes of mRNA metabolism such as pre-mRNA splicing. Its function in translational initiation regulation can either be enhanced by PAIP1 or repressed by PAIP2. PABPC1 can probably bind to cytoplasmic RNA sequences other than poly(A) in vivo. It may be involved in translationally coupled mRNA turnover. Implicated with other RNA-binding proteins in the cytoplasmic deadenylation/translational and decay interplay of the FOS mRNA mediated by the major coding-region determinant of instability (mCRD) domain.

产品性质

Alias: Polyadenylate-binding protein 1, PABP-1, Poly(A)-binding protein 1, PABPC1, PAB1, PABP1, PABPC2

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: P11940

Clone: -

Entrez GeneID: 26986

Host\Isotype: Rabbit IgG

Reactivity: Human, Mouse

Purification: affinity purified by Protein G

Applications: WB=1/1000;IHC=1/100-1/500

Immunogen: This PABPC1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 250-279 amino acids from the Central region of human PABPC1.

Formulation: Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol. prepared by Saturated Ammonium Sulfate (SAS) .

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.



网站: <https://www.aoruisai.com>

<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信